

Product datasheet for **TA399618L**

F3 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	F3 / Factor III / Tissue Factor / CD142
Formulation:	< 0.001 EU/μg, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144.66 kDa
Database Link:	P13726



[View online »](#)