

## Product datasheet for **TA399761L**

### **TMEFF1 Monoclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	TMEFF1 / Tomoregulin-1
Formulation:	< 0.001 EU/μg, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	150 kDa
Database Link:	<a href="#">Q8IYR6</a>



[View online »](#)