

Product datasheet for **TA399767**

TNFRSF4 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Macaque, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	TNFRSF4 / OX40 / CD134
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.22 kDa
Database Link:	P43489



[View online »](#)