

Product datasheet for **TA399786L**

FLT4 Monoclonal Antibody

Product data:

| | |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | Primary Antibodies |
| Applications: | Animal Model, ELISA, FACS, FN, Kinetics |
| Reactivity: | Human |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Specificity: | VEGFR3 / FLT4 |
| Formulation: | < 0.001 EU/μg, determined by LAL method. |
| Concentration: | 1 mg/ml |
| Purification: | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column. |
| Conjugation: | Unconjugated |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 146.18 kDa |
| Database Link: | P35916 |



[View online »](#)