

Product datasheet for **TA399873**

HAVCR2 Mouse Monoclonal Antibody [Clone ID: 7H2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	7H2
Applications:	ELISA, IHC, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	The original antibody was generated by immunizing mice with recombinant TIM-3 protein.
Specificity:	<p>The original antibody is specific for TIM-3. TIM-3 is a type I transmembrane protein and contains an immunoglobulin and a mucin-like domain in its extracellular portion and a tyrosine phosphorylation motif in its cytoplasmic portion.</p> <p>The antibody specifically detected TIM-3 protein in Western blot analysis. The antibody 7H2 and the commercial antibody D5D5R were used to stain tonsils and appendix tissue. Results showed the antibody specificity is comparable to the commercial antibody (CN113234158A).</p>
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Database Link:	Q8TDQ0
Synonyms:	CD366; CGMCCNO 20763; D5D5R; E3; HAVcr-2; Hepatitis A virus cellular receptor 2; T-cell immunoglobulin and mucin domain-containing protein 3; T-cell immunoglobulin mucin receptor 3; T-cell membrane protein 3; TIMD-3



[View online »](#)