

## Product datasheet for **TA399881**

### Mouse Monoclonal Antibody [Clone ID: 33D3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	33D3
Applications:	Dot, IF, MeDIP, MeDIP-seq
Reactivity:	Human, Mouse, Broad Range
Host:	Mouse
Clonality:	Monoclonal
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)