

Product datasheet for **TA399883**

H3C1 Escherichia coli Monoclonal Antibody

Product data:

| | |
|----------------|--|
| Product Type: | Primary Antibodies |
| Applications: | ChIP/ChIP-seq, IF |
| Reactivity: | Human, Mouse |
| Host: | Escherichia coli |
| Clonality: | Monoclonal |
| Concentration: | lot specific |
| Conjugation: | Biotin |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Database Link: | Q93081 |



[View online »](#)