

## Product datasheet for **TA399887**

### EP300 Mouse Monoclonal Antibody [Clone ID: 1H11G2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	1H11G2
Applications:	ChIP/ChIP-seq
Reactivity:	Human
Host:	Mouse
Clonality:	Monoclonal
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Database Link:	<a href="#">Q09472</a>



[View online »](#)