

Product datasheet for **TA420020**

Goat IgG (Fc) fragment Rabbit mAb Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:2000 - 1:6000
Reactivity:	Goat, Widely Range
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Goat IgG
Formulation:	Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	- Observed MW: 62kDa



[View online »](#)