

## Product datasheet for **TA420458**

### ADAMTS3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Propeptide domain
Formulation:	Aqueous solution sodium chloride (500 mM), sodium phosphate (10mM), Potassium Phosphate (1.8 mM) and Glycerol (50% v/v)
Concentration:	1.0 milligram/ml
Conjugation:	Unconjugated
Database Link:	<u><a href="#">Q15072</a></u>



[View online »](#)