

## Product datasheet for **TA420646**

### **MMP14 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Aminoterminal end of active enzyme
Formulation:	Aqueous solution sodium chloride (500 mM), sodium phosphate (10mM), Potassium Phosphate (1.8 mM) and Glycerol (50% v/v)
Concentration:	1.0 milligram/ml
Conjugation:	Unconjugated
Database Link:	<u><a href="#">P50281</a></u>



[View online »](#)