

Product datasheet for **TA421186**

VWF Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Carboxyterminal end of the ADAMTS-13 cleavage site
Formulation:	Aqueous solution sodium chloride (500 mM), sodium phosphate (10mM), Potassium Phosphate (1.8 mM) and Glycerol (50% v/v)
Concentration:	1.0 milligram/ml
Conjugation:	Unconjugated
Database Link:	<u>P04275</u>



[View online »](#)