

Product datasheet for **TA421315AM**

Anti-Mesothelin VHH-human IgG1 Fc, Biotinylated

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, FC
Recommended Dilution:	Flow Cyt 1:100
Reactivity:	Human
Host:	Human
Clonality:	Monoclonal
Formulation:	50% Glycerol; PBS, pH 7.4; 0.1% BSA Preservative: 0.1% Procline 300
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Concentration:	Lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Stability:	12 months from date of despatch
Database Link:	O00187
Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	MASP-2
Note:	Common Name: OMS721



Product images:




Figure 1. MASP2 protein is highly expressed on the surface of Jurkat cell membrane. Flow cytometry analysis with 15 µg/mL Anti-MASP2(narsoplimab biosimilar) mAb ([TA421315]) (Blue histogram) or isotype control mAb (Red histogram) on Jurkat cells.