

Product datasheet for **TA421349**

Anti-B7-H4(SGN-B7H4V biosimilar) mAb

Product data:

Product Type:	Primary Antibodies
Applications:	FC
Recommended Dilution:	Flow Cyt 1:100
Reactivity:	Human
Host:	Human
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Formulation:	50% Glycerol; PBS, pH 7.4; 0.1% BSA Preservative: 0.1% Procline 300
Reconstitution Method:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Concentration:	Lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.
Stability:	12 months from date of despatch
Database Link:	Q7Z7D3
Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	VTCN1
Note:	Common Name: SGNB7H4V,SGN-B7H4V,SGNB 7H4V




Product images:

Figure 1. Flow cytometry analysis with 1 µg/mL Anti-B7-H4(SGN-B7H4V biosimilar) mAb (TA421349) on HEK293 cells transfected with Human B7-H4 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).