

## Product datasheet for **TA421363AM**

### Anti-HEL Human IgG2-Kappa Isotype control mAb, Biotinylated

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, FC
Recommended Dilution:	ELISA 1:5000-10000 FC 1:100
Reactivity:	Ape
Host:	Chimeric
Isotype:	IgG2
Clonality:	Monoclonal
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.
Reconstitution Method:	1. Before opening the tube cap, centrifuge the sample tube at 5000g for 3-5min at room temperature to ensure the lyophilized sample settles down at the bottom of the tube. 2. Calculate the volume for reconstitution (in $\mu\text{L}$ ) using the formula: $[\text{Quantity (mg)}/\text{Concentration (mg/mL)}] \times 1000$ 3. Dissolve the lyophilized protein sample in sterile water based on the calculated volume ( $\mu\text{L}$ ) After adding sterile water, cover the lid and mix by gently tapping the tube 5-10 times. Note: Do not vortex or vigorously pipette the sample.
Concentration:	Lot specific
Conjugation:	Biotin
Stability:	12 months from date of dispatch
Database Link:	<u><a href="#">N/A</a></u>
Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.



[View online »](#)