

Product datasheet for TA421621


Entactin (NID1) Rabbit Monoclonal Antibody [Clone ID: R08-3H3]

Product data:


Product Type:	Primary Antibodies
Clone Name:	R08-3H3
Applications:	IHC, WB
Recommended Dilution:	WB: 1/500-1/1000 IHC: 1/50-1/100
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human Entactin/NID
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Predicted Protein Size:	Calculated MW: 136 kDa; Observed MW: 136 kDa
Gene Name:	nidogen 1
Database Link:	Entrez Gene 4811 Human P14543
Synonyms:	Enactin; Ent; Entactin; Entactin 1; NID; NID 1; Nid1; Nidogen; Nidogen 1

[View online »](#)

Product images:



Western blot analysis of Entactin/NID in HeLa, A549, HL-60, U251, U87-MG lysates using Entactin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Entactin/NID antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.