

Product datasheet for **TA421666**


CD53 Rabbit Monoclonal Antibody [Clone ID: R08-9P-5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	R08-9P-5
Applications:	FC, IHC, WB
Recommended Dilution:	WB:1/500-1000;IHC:1/50-100
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Peptide
Formulation:	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Predicted Protein Size:	Calculated MW:24 kDa;Observed MW: 35-55kDa
Gene Name:	CD53 molecule
Database Link:	Entrez Gene 963 Human P19397
Synonyms:	CD53; Cell surface glycoprotein CD53; Leukocyte surface antigen CD53; MOX44; Tetraspanin-25; Tspan-25; TSPAN25



Product images:



Western blot analysis of CD53 in Raji lysates using CD53 antibody.