

Product datasheet for **TA421721**

IL12A Rabbit Monoclonal Antibody [Clone ID: R08-1M-6]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R08-1M-6
Applications:	FC, IHC, WB
Recommended Dilution:	WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Peptide
Formulation:	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Predicted Protein Size:	Calculated MW:25 kDa;Observed MW: 35 kDa
Gene Name:	interleukin 12A
Database Link:	Entrez Gene 3592 Human P29459
Synonyms:	IL12A; Interleukin-12 subunit alpha, IL-12 subunit p35, IL-12A, Cytotoxic Lymphocyte Maturation Factor 35 kDa; NK cell stimulatory factor chain 1; NKSF1



[View online »](#)

Product images:



Western blot analysis of IL-12A in NIH3T3, Raw 264.7 lysates using IL-12A antibody.