

## Product datasheet for **TA421831**


### **D4 (ARHGDIB) Rabbit Monoclonal Antibody [Clone ID: R01-3K7]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	R01-3K7
Applications:	ICC/IF, IHC, WB
Recommended Dilution:	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human D4 GDI
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Predicted Protein Size:	Calculated MW: 23 kDa; Observed MW: 23 kDa
Gene Name:	Rho GDP dissociation inhibitor beta
Database Link:	<a href="#">Entrez Gene 397 Human P52566</a>
Synonyms:	D4; GDIA2; GDID4; Ly-GDI; LYGDI; RAP1GN1; RhoGDI2; Rho GDP dissociation inhibitor (GDI) beta

[View online »](#)

## Product images:



Western blot analysis of D4 GDI in Jurkat, U937 lysates using D4 GDI antibody.