

Product datasheet for **TA421848**


R Cadherin (CDH4) Rabbit Monoclonal Antibody [Clone ID: R04-5H6]

Product data:


| | |
|--------------------------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | R04-5H6 |
| Applications: | IP, WB |
| Recommended Dilution: | WB:1/500-1000 |
| Reactivity: | Human, Mouse, Rat |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Immunogen: | A synthetic peptide of human R Cadherin |
| Formulation: | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Concentration: | lot specific |
| Purification: | Affinity Purified |
| Conjugation: | Unconjugated |
| Storage: | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Predicted Protein Size: | Calculated MW: 100 kDa; Observed MW: 130 kDa |
| Gene Name: | cadherin 4 |
| Database Link: | Entrez Gene 1002 Human P55283 |
| Synonyms: | CAD4; cadherin 4; R-CAD; R-cadherin (retinal); RCAD; type 1 |



Product images:



Western blot analysis of RCadherin/Cadherin 4 in mouse brain, U87-MG lysates using R Cadherin antibody.



Western blot analysis of RCadherin/Cadherin 4 in rat Brain, C6 lysates using RCadherin/Cadherin 4 antibody.