

Product datasheet for **TA421852**

Gasdermin D (GSDMD) Rabbit Monoclonal Antibody [Clone ID: R07-2C4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	R07-2C4
Applications:	WB
Recommended Dilution:	WB: 1/500-1/1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human Cleaved N-terminal GSDMD
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Predicted Protein Size:	Calculated MW: 53 kDa
Gene Name:	gasdermin D
Database Link:	Entrez Gene 79792 Human P57764
Synonyms:	DFNA5L; FKSG10; Gasdermin-D; Gasdermin domain-containing protein 1; GSDMDC1

[View online »](#)