

Product datasheet for **TA422060**

Optineurin (OPTN) Rabbit Monoclonal Antibody [Clone ID: R01-2C5]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R01-2C5
Applications:	IHC, WB
Recommended Dilution:	WB: 1/500-1/1000 IHC: 1/50-1/100
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human Optineurin
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Predicted Protein Size:	Calculated MW: 66 kDa; Observed MW: 75 kDa
Gene Name:	optineurin
Database Link:	Entrez Gene 10133 Human Q96CV9
Synonyms:	ALS12; FIP2; GLC1E; HIP7; HYPL; NRP; TFIIA-INTP




[View online »](#)

Product images:



Western blot analysis of Optineurin in HeLa, A549, HL-60, U2OS, U87-MG lysates using Optineurin antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Optineurin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.