


Product datasheet for **TA422176**

MST3 (STK24) Rabbit Monoclonal Antibody [Clone ID: R03-6V-1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	R03-6V-1
Applications:	FC, ICC/IF, IHC, WB
Recommended Dilution:	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Peptide
Formulation:	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Predicted Protein Size:	Calculated MW:49 kDa;Observed MW: 49 kDa
Gene Name:	serine/threonine kinase 24
Database Link:	Entrez Gene 8428 Human Q9Y6E0
Synonyms:	HEL-S-95; MST3; MST3B; STE20; STK3


[View online »](#)

Product images:

Western blot analysis of MST3 in HeLa lysates using MST3 antibody.