


Product datasheet for TA423734

Cd52 Rabbit Monoclonal Antibody [Clone ID: OTIRA1553]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTIRA1553
Applications:	IP, WB
Recommended Dilution:	WB 1:500~1:2000 IP 1:50
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human CD52
Specificity:	CD52 Antibody detects endogenous levels of total CD52
Formulation:	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	Affinity-chromatography
Conjugation:	Unconjugated
Storage:	Store at +4°C short term. Store at -21°C long term. Avoid freeze / thaw cycle.
Predicted Protein Size:	8kDa
Gene Name:	CD52 antigen
Database Link:	Entrez Gene 23833 Mouse Q64389
Background:	May play a role in carrying and orienting carbohydrate, as well as having a more specific role.
Synonyms:	CAMPATH 1; CD52; CDw52; He5



Product images:

Western blot analysis of CD52 expression in
Mouse spleen cell lysate.