

Product datasheet for **TA424601S**

CAT Rabbit Monoclonal Antibody [Clone ID: 23GB2625]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB2625
Applications:	ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human Catalase
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	P04040
Synonyms:	CAT; Catalase; EC 1.11.1.6; Epididymis Secretory Sperm Binding Protein



[View online »](#)