

Product datasheet for **TA424794S**

GLA Rabbit Monoclonal Antibody [Clone ID: 23GB5545]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB5545
Applications:	FC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human Galactosidase alpha
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	P06280
Synonyms:	Agalsidase; Agalsidase Alfa; Alpha-D-Galactosidase A; Alpha-D-Galactoside Galactohydrolase; Alpha-D-Galactoside Galactohydrolase 1; Alpha-Gal A; Alpha-Galactosidase A; EC 3.2.1; EC 3.2.1.22; GALA; Galactosidase Alpha; Galactosylgalactosylglucosylceramidase GLA; GLA; Melibiase



[View online »](#)