

Product datasheet for **TA424825S**

GOLGA2 Rabbit Monoclonal Antibody [Clone ID: 23GB5570]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB5570
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human GM130
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q08379
Synonyms:	130 KDa Cis-Golgi Matrix Protein; DEDHMB; GM130; GM130 Autoantigen; GOLGA2; Golgi Autoantigen, Golgin Subfamily A, 2; Golgi Matrix Protein GM130; GOLGIN-95; Golgin-95; Golgin A2; Golgin Subfamily A Member 2; SY11 Protein



[View online »](#)