

Product datasheet for **TA424843S**

HDAC3 Rabbit Monoclonal Antibody [Clone ID: 23GB860]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB860
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human HDAC3
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q15379
Synonyms:	EC 3.5.1.-; EC 3.5.1.98; HD3 KDAC3; Histone Deacetylase 3; Protein Deacetylase HDAC3; Protein Deacylase HDAC3; RPD3; RPD3-2; SMAP45



[View online »](#)