

Product datasheet for **TA424880S**

IGFBP7 Rabbit Monoclonal Antibody [Clone ID: 24GB4910]

Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB4910
Applications:	FC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human IGFBP7
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q16270
Synonyms:	AGM; Angiomodulin; FSTL2; IBP-7; IGF-Binding Protein 7; IGFBP-7; IGFBP-7v; IGFBP-RP1; IGFBP7; IGFBPRP1; Insulin-Like Growth Factor-Binding Protein 7; Insulin-Like Growth Factor Binding Protein 7; Insulin Like Growth Factor Binding Protein 7; MAC25; MAC25 Protein; PGI2-Stimulating Factor; Prostacyclin-Stimulating Factor; PSF; RAMSVPS; TAF; Tumor-Derived Adhesion Factor



[View online »](#)