

Product datasheet for **TA424887S**

PPA1 Rabbit Monoclonal Antibody [Clone ID: 24GB1365]

Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB1365
Applications:	FC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human Pyrophosphatase 1
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q15181
Synonyms:	Cytosolic Inorganic Pyrophosphatase; Diphosphate Phosphohydrolase; EC 3.6.1.1; Epididymis Secretory Sperm Binding Protein Li 66p; HEL-S-66p; Inorganic Diphosphatase 1; Inorganic Pyrophosphatase; Inorganic Pyrophosphatase 1; PP; PP1; PPA1; PPase; Ppase; Pyrophosphatase (Inorganic); Pyrophosphatase (Inorganic)1; Pyrophosphatase 1; Pyrophosphate Phospho-Hydrolase; SID6-8061

[View online »](#)