

Product datasheet for **TA425112**

RPS6KA1 Rabbit Monoclonal Antibody [Clone ID: 23GB2425]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB2425
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human Phospho-RSK1(S380)
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q15418
Synonyms:	90 KDa Ribosomal Protein S6 Kinase 1; DJ590P13.1 (Ribosomal Protein S6 Kinase, 90kD, Polypeptide 1); EC 2.7.11; EC 2.7.11.1; HU-1; MAPK-Activated Protein Kinase 1a; MAPKAPK-1a; MAPKAPK1; MAPKAPK1A; MAPKAP Kinase 1a; MAP Kinase-Activated Protein Kinase 1a; P90-RSK 1; P90Rsk; P90RSK; P90RSK1; P90S6K; Ribosomal Protein S6 Kinase, 90kD, Polypeptide 1; Ribosomal Protein S6 Kinase, 90kDa, Polypeptide 1; Ribosomal Protein S6 Kinase A1; Ribosomal Protein S6 Kinase Alpha-1; Ribosomal S6 Kinase 1; RPS6KA1; RSK; RSK-1; RSK1; S6K-Alpha-1; S6K-Alpha 1



[View online »](#)