

## Product datasheet for **TA425158S**

### **PSMB10 Rabbit Monoclonal Antibody [Clone ID: 23GB2165]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	23GB2165
Applications:	FC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human PSMB10
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	<a href="#">P40306</a>
Synonyms:	BETA2i; Beta2i; EC 3.4.25.1; LMP10; Low Molecular Mass Protein 10; Macropain Subunit MECL-1; MECL1; MGC1665; Multicatalytic Endopeptidase Complex Subunit MECL-1; PRAAS5; Proteasome (Prosome, Macropain) Subunit, Beta Type, 10; Proteasome 20S Subunit Beta 10; Proteasome Catalytic Subunit 2i; Proteasome MECL-1; Proteasome Subunit Beta-2i; Proteasome Subunit Beta2i; Proteasome Subunit Beta 7i; Proteasome Subunit Beta 10; Proteasome Subunit Beta Type-10; Proteasome Subunit MECL1; Proteasome Subunit B2i; PSMB10



[View online »](#)