

Product datasheet for **TA425190S**

RBBP4 Rabbit Monoclonal Antibody [Clone ID: 24GB1240]

Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB1240
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human RbAp48
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q09028
Synonyms:	CAF-1 Subunit C; CAF-I 48 KDa Subunit; CAF-I P48; Chromatin Assembly Factor/CAF-1 P48 Subunit; Chromatin Assembly Factor 1 Subunit C; Chromatin Assembly Factor I P48 Subunit; Histone-Binding Protein RBBP4; Lin-53; MSI1 Protein Homolog; Nucleosome-Remodeling Factor Subunit RBAP48; NURF55; RBAP48; RbAp48; RB Binding Protein 4, Chromatin Remodeling Factor; RBBP-4; Retinoblastoma-Binding Protein 4; Retinoblastoma-Binding Protein P48; Retinoblastoma Binding Protein 4



[View online »](#)