

Product datasheet for **TA425221S**

SMAD5 Rabbit Monoclonal Antibody [Clone ID: 23GB2910]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB2910
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human Smad5
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q99717
Synonyms:	Dwfc; DWFC; HSmad5; JV5-1; MAD, Mothers Against Decapentaplegic Homolog 5; MAD, Mothers Against Decapentaplegic Homolog 5 (Drosophila); MADH5; MAD Homolog 5; Mothers Against Decapentaplegic, Drosophila, Homolog Of, 5; Mothers Against Decapentaplegic Homolog; Mothers Against DPP Homolog 5; SMA- And MAD-Related Protein 5; SMAD, Mothers Against DPP Homolog 5; SMAD, Mothers Against DPP Homolog 5 (Drosophila) 2; SMAD 5; Smad5; SMAD Family Member 5



[View online »](#)