

Product datasheet for **TA425273S**

ELOB Rabbit Monoclonal Antibody [Clone ID: 23GB5615]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB5615
Applications:	FC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human TCEB2
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q15370
Synonyms:	ELOB; EloB; Elongin, 18-KD Subunit; Elongin-B; Elongin 18 KDa Subunit; Elongin B; Epididymis Secretory Sperm Binding Protein; RNA Polymerase II Transcription Factor SIII P18 Subunit; RNA Polymerase II Transcription Factor SIII Subunit B; SIII; SIII P18; TCEB2; Transcription Elongation Factor B (SIII), Polypeptide 2 (18kDa, Elongin B); Transcription Elongation Factor B Polypeptide; Transcription Elongation Factor B Subunit 2



[View online »](#)