

Product datasheet for **TA425287**

TNFAIP8 Rabbit Monoclonal Antibody [Clone ID: 24GB2315]


Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB2315
Applications:	FC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human TNFAIP8
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	O95379
Synonyms:	GG2-1; Head And Neck Tumor And Metastasis-Related Protein; MDC-3.13; NDED; NF-Kappa-B-Inducible DED-Containing Protein; SCC-S2; SCCS2; TNF-Induced Protein GG2-1; TNFAIP8; TNF Alpha-Induced Protein 8; TNF Alpha Induced Protein 8; Tumor Necrosis Factor, Alpha-Induced Protein 8; Tumor Necrosis Factor, Alpha Induced Protein 8; Tumor Necrosis Factor Alpha-Induced Protein 8




[View online »](#)


Product images:



Western blotting analysis using anti-TNF alpha induced protein 8 antibody . Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TNF alpha induced protein 8 antibody and HRP-conjugated goat anti-rabbit secondary antibody respectively. Image was developed using anti-FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-TNF alpha induced protein 8 antibody . TNF alpha induced protein 8 expression in wild type (WT) and TNF alpha induced protein 8 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates . β-Tubulin serves as a loading control. The blot was incubated with anti-TNF alpha induced protein 8 antibody and HRP-conjugated goat anti-rabbit secondary antibody respectively. Image was developed using anti-NaQ™ ECL Substrate Kit .



Flow cytometric analysis of TNF alpha induced protein 8 expression in Hela cells using anti-TNF alpha induced protein 8 antibody . Green, isotype control; red, TNF alpha induced protein 8.