

Product datasheet for **TA425293S**

TCF7L1 Rabbit Monoclonal Antibody [Clone ID: 24GB1985]

Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB1985
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human TCF3
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q9HCS4
Synonyms:	HMG Box Transcription Factor 3; TCF-3; TCF3; TCF7L1; Transcription Factor 7 Like 1; Transcription Factor 7-Like 1; Transcription Factor 7-Like 1 (T-Cell Specific, HMG-Box)



[View online »](#)