

## Product datasheet for **TA425311S**

### UBE2T Rabbit Monoclonal Antibody [Clone ID: 24GB2350]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB2350
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human HSPC150
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	<a href="#">Q9NPD8</a>
Synonyms:	Cell Proliferation-Inducing Gene 50 Protein; E2 Ubiquitin-Conjugating Enzyme T; EC 2.3.2.23; FANCT; HSPC150; HSPC150 Protein Similar To Ubiquitin-Conjugating Enzyme; PIG50; UBE2T; Ubiquitin-Conjugating Enzyme E2 T; Ubiquitin-Conjugating Enzyme E2T; Ubiquitin-Conjugating Enzyme E2T (Putative); Ubiquitin-Protein Ligase T; Ubiquitin Carrier Protein T; Ubiquitin Conjugating Enzyme E2 T; Ubiquitin Conjugating Enzyme E2T



[View online »](#)