

Product datasheet for **TA425323S**

UFC1 Rabbit Monoclonal Antibody [Clone ID: 23GB4245]

Product data:

Product Type:	Primary Antibodies
Clone Name:	23GB4245
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human UFC1
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q9Y3C8
Synonyms:	HSPC155; NEDSG; Ubiquitin-Fold Modifier-Conjugating Enzyme 1; Ubiquitin-Fold Modifier Conjugating Enzyme 1; UFC1; Ufm1-Conjugating Enzyme 1



[View online »](#)