

Product datasheet for **TA425344S**

ZPR1 Rabbit Monoclonal Antibody [Clone ID: 24GB1805]

Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB1805
Applications:	FC, ICC, WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000, Flow Cytometry (FCM): 1:2,000, Immunocytochemistry (IC): 1:100-1:1,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthesized peptide derived from human ZNF259
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	Q75312
Synonyms:	GKAF; Zinc Finger Protein 259; Zinc Finger Protein ZPR1; ZNF259; ZPR1; ZPR1 Zinc Finger



[View online »](#)