

Product datasheet for **TA425349S**

AK1 Mouse Monoclonal Antibody [Clone ID: 24GB2940]

Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB2940
Applications:	WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human Adenylate Kinase 1
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	P00568
Synonyms:	Adenylate Kinase 1; Adenylate Kinase Isoenzyme 1; Adenylate Monophosphate Kinase; AK 1; AK1; ATP-AMP Transphosphorylase 1; ATP:AMP Phosphotransferase; EC 2.7.4; EC 2.7.4.3; EC 2.7.4.6; EC 2.7.4.10; Epididymis Secretory Sperm Binding Protein; HTL-S-58j; Myokinase; Testis Secretory Sperm Binding Protein Li 58j

[View online »](#)