

Product datasheet for **TA425366S**

ANXA1 Mouse Monoclonal Antibody [Clone ID: 24GB6985]

Product data:

Product Type:	Primary Antibodies
Clone Name:	24GB6985
Applications:	WB
Recommended Dilution:	Western Blotting (WB): 1:1,000-1:5,000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human Annexin A1 beta
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at -20 °C for one year.
Database Link:	P04083
Synonyms:	Annexin-1; Annexin A1; Annexin I; Annexin I (Lipocortin I); ANX1; ANXA1; Calpactin-2; Calpactin II; Chromobindin-9; Epididymis Secretory Sperm Binding Protein; Lipocortin I; LPC1; P35; Phospholipase A2 Inhibitory Protein



[View online »](#)