

Product datasheet for **TA425374**


RAB27A Mouse Monoclonal Antibody [Clone ID: 24GB7165]

Product data:



Product Type:	Primary Antibodies
Clone Name:	24GB7165
Applications:	WB
Recommended Dilution:	Western Blotting (WB): 1:1,000–1:5,000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human RAB27
Formulation:	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Concentration:	Lot dependent
Purification:	Affinity Purified
Conjugation:	Unconjugated
Stability:	Store at –20 °C for one year.
Database Link:	P51159
Synonyms:	EC 3.6.5.2; GS2; GTP-Binding Protein Ram; HsT18676; Mutant Ras-Related Protein Rab-27A; Rab-27; RAB27; RAB27A; RAB27A, Member RAS Oncogene Family; RAM; Ras-Related Protein Rab-27A



Product images:



Western blotting analysis using anti-RAB27A antibody . Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-RAB27A antibody and HRP-conjugated goat anti-mouse secondary antibody respectively. Image was developed using anti-FeQ™ ECL Substrate Kit .



Western blotting analysis using anti-RAB27A antibody . RAB27A expression in wild-type (WT) and RAB27A shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-RAB27A antibody and HRP-conjugated goat anti-mouse secondary antibody respectively. Image was developed using anti-FeQ™ ECL Substrate Kit .