

## Product datasheet for **TA500001**

### **AKT3 Mouse Monoclonal Antibody [Clone ID: OTI9B2]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI9B2
<b>Applications:</b>	FC, IF, IHC, WB
<b>Recommended Dilution:</b>	WB 1:500, IHC 1:50, IF 1:100, FLOW 1:100
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG3
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human AKT3 (NP_005456) produced in E.coli.
<b>Formulation:</b>	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	0.25 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Predicted Protein Size:</b>	55.6 kDa
<b>Gene Name:</b>	AKT serine/threonine kinase 3
<b>Database Link:</b>	<a href="#">NP_005456</a> <a href="#">Entrez Gene 23797 Mouse</a> <a href="#">Entrez Gene 29414 Rat</a> <a href="#">Entrez Gene 10000 Human</a> <a href="#">Q9Y243</a>
<b>Background:</b>	The protein encoded by this gene is a member of the AKT, also called PKB, serine/threonine protein kinase family. AKT kinases are known to be regulators of cell signaling in response to insulin and growth factors. They are involved in a wide variety of biological processes including cell proliferation, differentiation, apoptosis, tumorigenesis, as well as glycogen synthesis and glucose uptake. This kinase has been shown to be stimulated by platelet-derived growth factor (PDGF), insulin, and insulin-like growth factor 1 (IGF1). Alternatively splice transcript variants encoding distinct isoforms have been described.




[View online »](#)

**Synonyms:** MPPH; MPPH2; PKB-GAMMA; PKBG; PRKBG; RAC-gamma; RAC-PK-gamma; STK-2


**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase

**Protein Pathways:** Acute myeloid leukemia, Adipocytokine signaling pathway, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, MAPK signaling pathway, Melanoma, mTOR signaling pathway, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Tight junction, Toll-like receptor signaling pathway, VEGF signaling pathway


**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY AKT3 (Cat# [RC221051], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AKT3 (Cat# TA500001).



AKT3 antibody ( 9B2 ) at 1:500 dilution + BALB/C 3T3 cell lysate



Western blot analysis of extracts (10ug) from 2 different cell lines by using anti-AKT3 monoclonal antibody at 1:200 dilution.






Figure from citation: Western Blot of Akt3 protein level by using anti-Akt3 antibody in rat endometrium during pregnancy and pseudopregnancy. [View Citation](#)



Immunohistochemical staining of paraffin-embedded Human lung tissue within the normal limits using anti-AKT3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500001)



Anti-AKT3 mouse monoclonal antibody (TA500001) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY AKT3 ([RC221051]).



HEK293T cells transfected with either overexpress plasmid ([RC221051], Red) or empty vector control plasmid (Blue) were immunostained by anti-AKT3 antibody (TA500001), and then analyzed by flow cytometry.