

Product datasheet for TA500021AM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Cytokeratin 8 (KRT8) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1B12]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1B12
Applications: IF, IHC, WB

Recommended Dilution: WB 1:2500~5000, IHC 1:50 - 1:150, IF 1:50 - 1:100

Reactivity: Human, Dog, Rat, Monkey, Mouse

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 91-381 of human CK8

(NP_002264) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 53.5 kDa

Gene Name: keratin 8

Database Link: NP 002264

Entrez Gene 16691 MouseEntrez Gene 25626 RatEntrez Gene 486513 DogEntrez Gene 3856

<u>Human</u> P05787



Cytokeratin 8 (KRT8) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1B12] – TA500021AM


Background:

Keratin 8 is a member of the type II keratin family clustered on the long arm of chromosome 12. Type I and type II keratins heteropolymerize to form intermediate-sized filaments in the cytoplasm of epithelial cells. The product of this gene typically dimerizes with keratin 18 to form an intermediate filament in simple single-layered epithelial cells. This protein plays a role in maintaining cellular structural integrity and also functions in signal transduction and cellular differentiation. Mutations in this gene cause cryptogenic cirrhosis.


Synonyms: CARD2; CK-8; CK8; CYK8; K2C8; K8; KO

Protein Families: Druggable Genome

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY KRT8 (Cat# [RC209570], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-KRT8 (Cat# [TA500021])(1:500).




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-KRT8 monoclonal antibody.






Western blot analysis of extracts (10ug) from 8 Human tissue by using anti-KRT8 monoclonal antibody at 1:1000 (1: Testis; 2: Uterus; 3: Breast; 4: Brain; 5: Liver; 6: Ovary; 7: Thyroid gland; 8: Colon).




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-KRT8 mouse monoclonal antibody ([TA500021]) at 1:150 dilution.



Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-KRT8 mouse monoclonal antibody ([TA500021]) at 1:150 dilution.






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-KRT8 mouse monoclonal antibody ([TA500021]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-KRT8 mouse monoclonal antibody ([TA500021]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-KRT8 mouse monoclonal antibody ([TA500021]) at 1:150 dilution.


Cytokeratin 8 (KRT8) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1B12] – TA500021AM



Immunohistochemical staining of paraffinembedded Human bladder tissue within the normal limits using anti-KRT8 mouse monoclonal antibody ([TA500021]) at 1:150 dilution.




Anti-KRT8 mouse monoclonal antibody ([TA500021]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY KRT8 ([RC209570]).



Immunofluorescent staining of HeLa cells using anti-KRT8 mouse monoclonal antibody ([TA500021]).





Immunofluorescent staining of A549 cells using anti-KRT8 mouse monoclonal antibody ([TA500021]).