

Product datasheet for TA500030AM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Neurotrophin 3 (NTF3) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5A2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5A2
Applications: IHC, WB

Recommended Dilution: WB: 1:5000, IHC 1:50, IF (1:100)

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Recombinant protein expressed in E.coli corresponding to amino acids 139-257 of human

NT3

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 27.4 kDa

Gene Name: neurotrophin 3

Database Link: NP 001096124

Entrez Gene 18205 MouseEntrez Gene 81737 RatEntrez Gene 4908 Human

P20783



Neurotrophin 3 (NTF3) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5A2] – TA500030AM

Background: NT3 is a member of the neurotrophin family, that controls survival and differentiation of


mammalian neurons. It is closely related to both nerve growth factor and brain-derived neurotrophic factor. It may be involved in the maintenance of the adult nervous system, and may affect development of neurons in the embryo when it is expressed in human placenta. NT3-deficient mice generated by gene targeting display severe movement defects of the limbs. The mature peptide of this protein is identical in all mammals examined including human, pig, rat and mouse

Synonyms: HDNF; NGF-2; NGF2; NT-3; NT3


Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: MAPK signaling pathway, Neurotrophin signaling pathway

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NTF3 ([RC214281], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NTF3. Positive lysates [LY420181] (100ug) and [LC420181] (20ug) can be purchased separately from OriGene.




NT3 antibody (5A2) at 1:5000 + Mature Recombinant Human NT3 (14KDa) (Cat# [RC214281])


Neurotrophin 3 (NTF3) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5A2] – TA500030AM




Immunohistochemical staining of paraffinembedded breast within the normal limits using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).




Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).




Immunohistochemical staining of paraffinembedded kidney within the normal limits using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).



Immunohistochemical staining of paraffinembedded Carcinoma of kidney using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).




Immunohistochemical staining of paraffinembedded lung within the normal limits using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).




Immunohistochemical staining of paraffinembedded Carcinoma of lung using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).

Neurotrophin 3 (NTF3) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5A2] – TA500030AM



Immunohistochemical staining of paraffinembedded Carcinoma of pancreas using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min).



Immunohistochemical staining of paraffinembedded endometrium within the normal limits using anti-NTF3 ([TA500030]) mouse monoclonal antibody (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min)