

### **Product datasheet for TA500212S**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Cytokeratin 19 (KRT19) Mouse Monoclonal Antibody [Clone ID: OTI3F8]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI3F8

**Applications:** IF, IHC, WB

**Recommended Dilution:** WB 1:1000, IHC 1:50 - 1:150, IF 1:100

Reactivity: Human, Dog, Rat, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Recombinant fragment expressed in E.coli corresponding to amino acids 240-390 of human

CK19

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.15 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 43.9 kDa

Gene Name: keratin 19

Database Link: NP 002267

Entrez Gene 360626 RatEntrez Gene 100685481 DogEntrez Gene 3880 Human

P08727







### Background:


Keratin 19 is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. Unlike its related family members, this smallest known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically expressed in the periderm, the transiently superficial layer that envelopes the developing epidermis.

**Synonyms:** CK19; K1CS; K19

# **Product images:**









Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-KRT19 monoclonal antibody.







Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-KRT19 mouse monoclonal antibody ([TA500212]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-KRT19mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500212], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Carcinoma of kidney tissue using anti-KRT19mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500212], Dilution 1:50)




Immunohistochemical staining of paraffinembedded liver tissue within the normal limits using anti-KRT19mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500212], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-KRT19 mouse monoclonal antibody ([TA500212]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of endometrium tissue using anti-KRT19mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500212], Dilution 1:50)



Immunohistochemical staining of paraffinembedded Human bladder tissue within the normal limits using anti-KRT19 mouse monoclonal antibody ([TA500212]) at 1:150 dilution.



Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-KRT19 mouse monoclonal antibody ([TA500212]) at 1:150 dilution.



Anti-KRT19 mouse monoclonal antibody ([TA500212]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY KRT19 ([RC209707]).