

Product datasheet for TA500274S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Flt3 ligand (FLT3LG) Mouse Monoclonal Antibody [Clone ID: OTI8D10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI8D10
Applications: IF, WB

Recommended Dilution: WB 1:1000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Recombinant protein expressed in E.coli corresponding to amino acids 27-181 of human

FLT3LG

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 26.4 kDa

Gene Name: fms related receptor tyrosine kinase 3 ligand

Database Link: NP 001450

Entrez Gene 2323 Human

<u>P49771</u>

Background: Stimulates the proliferation of early hematopoietic cells. Synergizes well with a number of

other colony stimulating factors and interleukins

Synonyms: FL; FLT3L


Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein, Transmembrane


Protein Pathways: Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Pathways in cancer





Product images:





HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY FLT3LG (Cat# [RC222242], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FLT3LG(Cat# [TA500274]). Positive lysates [LY419927] (100ug) and [LC419927] (20ug) can be purchased separately from OriGene.

Anti-FLT3LG mouse monoclonal antibody ([TA500274]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FLT3LG ([RC222242]).