

# **Product datasheet for TA500366S**

# OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### **SNAIL (SNAI1) Mouse Monoclonal Antibody [Clone ID: OTI10D7]**

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI10D7

Applications: WB

Recommended Dilution: WB 1:1000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human SNAI1 expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Conjugation:** Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 29.1 kDa

**Gene Name:** snail family transcriptional repressor 1

Database Link: NP 005976

Entrez Gene 20613 MouseEntrez Gene 116490 RatEntrez Gene 6615 Human

095863

**Background:** The Drosophila embryonic protein snail is a zinc finger transcriptional repressor which

downregulates the expression of ectodermal genes within the mesoderm. The nuclear

protein encoded by this gene is structurally similar to the Drosophila snail protein, and is also

thought to be critical for mesoderm formation in the developing embryo. At least two

variants of a similar processed pseudogene have been found on chromosome 2

Synonyms: dj710H13.1; SLUGH2; SNA; SNAH; SNAIL; SNAIL1







**Protein Families:** Druggable Genome

**Protein Pathways:** Adherens junction

# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY SNAI1 (Cat# [RC204581], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SNAI1(Cat# [TA500366]). Positive lysates [LY401811] (100ug) and [LC401811] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (10ug) from 1 cell line by using anti-SNAI1 monoclonal antibody at 1:200 dilution.