

## Product datasheet for TA500381BM

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### TRPM4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI10H5]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI10H5
Applications: IHC, WB

Recommended Dilution: WB 1:5000~10000, IHC 1:50

Reactivity: Human, Dog, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human TRPM4 (NP\_060106) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 134.3 kDa

Gene Name: transient receptor potential cation channel subfamily M member 4

Database Link: NP 060106

Entrez Gene 484385 DogEntrez Gene 719120 MonkeyEntrez Gene 54795 Human

O8TD43






Background:

Calcium-activated non selective (CAN) cation channel that mediates membrane depolarization. While it is activated by increase in intracellular Ca(2+), it is impermeable to it. Mediates transport of monovalent cations (Na(+) > K(+) > Cs(+) > Li(+)), leading to depolarize the membrane. It thereby plays a central role in cadiomyocytes, neurons from entorhinal cortex, dorsal root and vomeronasal neurons, endocrine pancreas cells, kidney epithelial cells, cochlea hair cells etc. Participates in T-cell activation by modulating Ca(2+) oscillations after T lymphocyte activation, which is required for NFAT-dependent IL2 production. Involved in myogenic constriction of cerebral arteries. Controls insulin secretion in pancreatic betacells. May also be involved in pacemaking or could cause irregular electrical activity under conditions of Ca(2+) overload.


Synonyms: hTRPM4; LTrpC4; PFHB1B; TRPM4B

**Protein Families:** Druggable Genome, Ion Channels: Transient receptor potential, Transmembrane

# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY TRPM4 (Cat# [RC216888], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRPM4(Cat# [TA500381]). Positive lysates [LY413636] (100ug) and [LC413636] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-TRPM4 monoclonal antibody.








Figure from citation: Western Blot of TRPM4 protein level by using anti-TRPM4 antibody in human endothelial cells. <u>View Citation</u>




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-TRPM4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500381])




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-TRPM4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500381])






Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-TRPM4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500381])



Immunohistochemical staining of paraffinembedded Human Ovary tissue within the normal limits using anti-TRPM4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500381])



Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-TRPM4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500381])